

Salinity Refractometer, 0 to 28% with ATC Model R9600



TECHNICAL DATA





Features

- Measures the concentration of dissolved salt in water (O to 28% Salinity)
- · Easy to operate with simple to read scales
- Requires only 2 or 3 drops of sample
- Automatic Temperature Compensation (ATC)
- Prism and lens with a simple focus adjustment provides repeatable results
- Includes cover plate, calibration screwdriver, cleaning cloth, solution/sample applicator and hard carrying case

Specifications

Measuring Range: O to 28% (Salinity) Accuracy: ±0.2% 0.20% Resolution:

Automatic Temperature Compensation Range:

Dimensions:

Weight:

Model

R9600

10 to 30°C 6.5 x 1 x 1.5"

(170 x 30 x 40 mm)

Description	
Salinity Refractometer, O to 28% with ATC	

7oz (200g)

RPDPA1 Replacement Refractometer Lens Cover Salinity Refractometer, O to 28% with ATC & NIST **R9600-NIST**